### Vescent FFC Frequency Combs

The FFC-100 is a fully stabilizable, octave-spanning frequency comb with precise control over  $f_{\rm rep}$ ,  $f_{\rm opt}$ , and  $f_{\rm CEO}$ . The Er-doped fiber MOPA architecture is simple and robust, and yet high-performance. A highly non-linear fiber broadens the spectrum and our unique  $f_{\rm CEO}$  lock detection reduces the size, weight, and power of the system. The FFC frequency combs are designed and built to ensure stable, low-phase noise operation with Allan Deviations supporting the next generation of optical atomic clocks, gravimeters, quantum computers, optical sensors, and more.

The FFC-100 was designed for low SWaP and turn-key, stable operation: A single 2U 19" rack mount chassis contains the oscillator, amplifier, pump lasers, supercontinuum generation module, and  $f_{\text{CEO}}$  detection and



FFC-100 turnkey frequency combs

lock as well as the control electronics. The simple oscillator mode locks at start up every time and the innovative passive SESAM mode-locker is specially designed for a robust, long life. Our unique oscillator design also makes it easy to precisely factory match the repetition rate of two (or more) FFC-100 combs for multi-comb spectroscopy experiments.

#### Features:

- Turn-key operation
- 1560 nm center wavelength
- Low phase & amplitude noise
- 2U 19" rack-mounted enclosure or modular
- ullet  $f_{\text{rep}}$  monitoring, control, and matching
- Input port for  $f_{\text{opt}}$
- Repetition rates from 80 to 250 MHz
- Optional visible extensions
- Optional super continuum flattening
- Made in America

### <u>Applications:</u>

- Time & frequency readout & transfer
- Frequency ruler
- Dual- and multi-comb spectroscopy
- Quantum sensing, computing, & cryptography
- Low-phase noise rf generation




Also available in compact modular form for OEM integration.




## Meet the FFC-100


Touch screen for easy set up, control, and monitoring Rep rate monitor Oscillator optical Supercontinuum monitor optical output Mode Lock status  $f_{\rm CEO}$  monitor indicator (3.5") $f_{
m opt}$  monitor Amplifier optical Rotary knob allows fast, precise Optical reference monitor value entry input 19" fcso Modified Allan Deviation 0



Full octave-spanning spectrum with optional flattened super-continuum. Dashed line indicates  $f_{\rm opt}$  minimum lockable power.



Rock-solid performance of the FFC offers favorable stability with respect to the next generation of atomic clock requirements.



Non-linear extension of comb teeth allows for referencing of  $f_{\rm opt}$  in the visible. 700-740 nm also available, but not shown for clarity.





# FFC-100 Specifications

Parameter	Units	Min.	Тур.	Max.	Comments
Center Wavelength	nm	1555	1565	1575	
Repetition Rate	MHz	100		200	Standard: 100,125,200 MHz
Optical Powers					All connectors PM FC/APC
Oscillator Average Output Power	mW	0.01			>15 nm optical bandwidth
Amplifier Average Output Power	mW	5			>30 nm optical bandwidth
Supercontinuum Output Power	mW	20			1000-2000 nm spectrum
RF Outputs					All connectors SMA
$f_{ extsf{CEO}}$ Signal-to-Noise Ratio	dВ	35			100 kHz RBW
$f_{ extsf{CEO}}$ Integrated Phase Noise	rad			1	10 Hz-1 MHz
$f_{ extsf{CEO}}$ Frequency Stability $^{^2}$	$ au^{ ext{-}1}$			5x10 <sup>-17</sup>	In-loop Allan Deviation
$f_{ m opt}$ Optical Input Power	mW			1	
$f_{ m opt}$ Signal-to-Noise Ratio $^{^3}$	dВ		40		100 kHz RBW
$f_{ m opt}$ Integrated Phase Noise $^3$	mrad		300		10 Hz-1 MHz
$f_{ m opt}$ Frequency Stability $^{2,3}$	$ au^{ ext{-}1}$		2.5x10 <sup>-17</sup>		In-loop Allan Deviation
$f_{rep}$ Output Power Level	dBm		0		
Frequency Transducers					
$f_{ extsf{CEO}}$ Tuning Range	MHz	100			Pump current tuning
$f_{ extsf{CEO}}$ Input Voltage Tuning Range	V	-10			SMA input
$f_{ m rep}$ PZT Tuning Range	Hz	30			Depends on repetition rate
$f_{ m rep}$ Input Voltage Tuning Range	V	0			SMA input
$f_{rep}$ Temperature Tuning Range	kHz	25			Touchscreen control 4
$f_{rep}$ Temperature Tuning Sensitivity	kHz/°C	1			
Size, Weight, and Power					
Line Voltage	VAC	90		250	50/60 Hz
Power Consumption	W			50	
Weight	kg		5.5		
Size	L		12.3		2U rack mount, 43.5 cm Width x 10.5 cm Height x 27 cm Depth

 $<sup>^{1}\</sup>mbox{Does}$  not include the SLICE-FPGA controller, which is sold separately.

<sup>&</sup>lt;sup>4</sup> Allows access to a slow servo correction mode that automatically changes oscillator cavity temperature so that the PZT is centered in its operating range.



 $<sup>^2\,\</sup>mbox{Allan}$  Deviation from zero-dead-time L counter with 1s gate time.

<sup>&</sup>lt;sup>3</sup>Depends on user-supplied optical reference. Data given for 1 kHz 1560 nm reference laser with > 0.1 mW input power when phase locked with a SLICE-FPGA.